PRODUCTION INFORMATION



PHOSPHATE ELISA DILUENT/BLOCKING BUFFER (10X)

Catalog No.: RE-DILU-P-100

Description: BSA in Phosphate buffered saline (10X)

Buffer: 10 mM phosphate, 150 mM NaCl, 0.05% Tween-20, 1 % BSA at pH 7.8-

8.0 (after dilution).

Size: 100 mL

Application: This assay diluent can be used in peroxidase-based ELISA or blotting.

Storage: Short term store at 2-8 °C (up to 12 months).

Instruction: Make 10-fold dilution of the 10X assay diluent buffer (1 part buffer with 9

part of distilled or de-ionized water). Working solution (after dilution) can

be stored at room temperature for up to 2 weeks.